




SHENZHEN XEAST TECHNOLOGY CO., LTD.

◆ We specialize in R&D, production and sales of various Laser Level, Thermal Imager, Laser Distance Meter, and other environmental equipment.
 ◆ We do our best to provide sincere pre-sales, on-sale and after-sales service to build a bright future together with you.

 +86-(0)755-82566735
  www.xeast.cn www.laser-level-supplier.com
  info@xeast.cn


 Quality assurance


 Reasonable price


 Fast delivery


 7/24 Online service

● SPECIFICATION



Medical supplies baby Infrared Digital Body Non-contact IR Infrared Thermometer

1. Introduction

This is a professional hand-held non-contact infrared thermometer which can reduce cross infection effectively. It is widely used for measuring temperature of human body and animals. Use the meter to test the temperature of temporal artery, which is the middle of the forehead from the hair line down to brows, could get stable and accurate measuring. Too much to right or left in the middle of the forehead will bring an inaccurate temperature reading.

2. Features

- * Precise non-contact measurement
- * Select-able temperature units °C/°F
- * Automatic measurement range selection with resolution of 0.1°C/ 0.1°F
- * LCD backlight
- * Auto data hold
- * Auto power off

3. Safety

- * Please read this user's manual thoroughly before using this unit
- * For the most accurate results, making reading with and ambient temperature of 10~40°C
- * Do not use the meter near electrified objects to avoid electric shock
- * Do not use the meter in the relative humidity over 85%

* Do not use volatile liquid to clean the unit, swipe it with dry soft cloth.

* Do not disassemble the unit, repair it by qualified personnel

4. Operating Instruction

4.1 Notice before operation

* If the meter has been stored in a cold or hot environment, allow it 15 to 20 minutes to acclimate to room temperature before making measurements.

* Measurement readings are varied from skin color and thickness

* Do not keep the meter in the sun or immerse it in water

* Make sure the skin is dry and that no hair interferes with the measurement.

Measuring parts	Normal Temperature
Anus	36.6°C~38.0°C/97.8°F~100.4°F
Oral	35.5°C~37.5°C/95.9°F~99.5°F
Armpit	34.7°C~37.3°C/94.4°F~99.1°F
Ear	35.8°C~38.0°C/96.4°F~100.4°F

4.2 Normal body temperature range for different parts of body

Ages	Normal Temperature
0~2	36.4°C~38.0°C/97.5°F~100.4°F
3~10	36.1°C~37.8°C/97.0°F~100.0°F
11~65	35.9°C~37.6°C/96.6°F~99.7°F

5. Specifications

Measuring Range	32.0°C~42.0°C(89.6°F~107.6°F)
Accuracy	32.0°C~35.9°C(89.6°F~96.6°F);±0.3°C/0.5°F 36.0°C~39.0°C(96.8°F~102.2°F);±0.2°C/0.4°F 39.1°C~42.0°C(102.3°F~107.6°F);±0.3°C/0.5°F
Resolution	0.1°C/0.1°F
Measuring Distance	5cm
Operating Temp.	10°C~50°C
Relative Humidity	≤85%
Auto Power Off	Meter shuts off automatically after 7 seconds of inactivity
Power Supply	9V battery
Weight	145g
Dimensions(L*W*H)	134*88.5*36mm
Note: the accuracy is obtained under ambient temperature 18°C-28°C/64°F-82°F.	





● APPLICATIONS



A smiling woman with long brown hair is sitting on a light-colored couch, holding a baby on her lap. The baby is wearing a striped shirt and is looking towards the camera. The background is a bright, airy room with a window and some furniture.

Take Your Own Temperature At Home Without Going Out

Provides simple, safe, and rapid body
temperature testing for adult health
and infant health



● **PACKAGE INCLUDE**



Package



Body thermometer



Color box



Battery



Manual



Shenzhen XEAST Technology Co., Ltd (XEAST TECHNOLOGY CO., LIMITED) is a leading supplier of Test & Measurement Instruments in China, such as Laser Level, Thermal Imaging Camera, Laser Distance Meter, Air Quality Tester, Sound level Meter, Infrared Thermometer, Anemometer, Digital Multimeter, Digital Clamp Meter and so on.

2015, our own brand-- XEAST kick off, our factory passed certificate of ISO9000:2000 Quality Control System, and provides CE, FCC and POHS-marked with XEAST brand to customers outside the mainland and in over 100 countries across America, Europe and Asia. With the advantage of high-grade products quality, competitive price and marvelous after-sales service, Laesent international enjoys great reputation at home and abroad.

The concept of Customer first, serving foremost Customers is the most important asset of our company, continuously trying to improve the products and service and devoting to responding to customers appeals and satisfaction timely so as to keep our customers loyalty and attract more potential customers. New century, new journey, XEAST is sure to keep pace with the advance of the times and have a free mind, build a strong teamwork to compete the challenges of new global times. XEAST sincerely expects to establish a goodwill and business relation with business people all of the world on the basis of mutual equality and benefit. We welcome all inquiries, whether for small or large transactions, knowing that step by step, together, we could build a better tomorrow!




CERTIFICATE

HONOR RECOGNITION IS TRUSTWORTHY



PATENTS



PAYMENT METHOD



PACKAGING & SHIPPING

1. DHL / UPS / FedEx / TNT , Door-to-Door.
2. By Air or by Sea , for FCL; Airport/ Port receiving.
3. Customers Specifying Freight Forwarders or Negotiable Shipping Methods.
4. We choose the best and safe packaging material to make sure your orders won't be damaged during delivery.





Company Outlook



Sales Dept



Showroom



Company View



Manufacturing Line



Warehouse



Testing Machine



SMT Dept



Mold Insection Line



Team Building



Customer Visting



